

Antibodies Product Datasheet

Product Name:	EXTL1 Antibody
Catalog Number:	C42189
Reactivity:	Human, Mouse
Applications:	WB
Dilutions:	WB 1:1000-2000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human EXTL1
Concentration:	1 mg/ml
Gene Name:	EXTL1 EXTL
Protein Name:	Exostosin-like 1 (EC 2.4.1.224) (Exostosin-L) (Glucuronosyl-N-acetylglucosaminyl-proteoglycan 4-alpha-N-acetylglucosaminyltransferase) (Multiple exostosis-like protein)
Other Name:	
Human Swiss Prot No.:	Q92935
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of human EXTL1
Storage and Stability:	-20°C/1 year