

Antibodies Product Datasheet

Product Name:	FBW10 Antibody
Catalog Number:	C42199
Reactivity:	Human, Mouse
Applications:	WB
Dilutions:	WB 1:1000-2000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human FBW10
Concentration:	1 mg/ml
Gene Name:	FBXW10
Protein Name:	F-box/WD repeat-containing protein 10 (F-box and WD-40 domain-containing protein 10) (Ubiquitin ligase-specificity factor)
Other Name:	
Human Swiss Prot No.:	Q5XX13
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of human FBW10
Storage and Stability:	-20°C/1 year