

Antibodies Product Datasheet

Product Name:	FIS1 Antibody
Catalog Number:	C42202
Reactivity:	Human, Mouse, Rat
Applications:	WB
Dilutions:	WB 1:1000-2000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human FIS1
Concentration:	1 mg/ml
Gene Name:	FIS1 TTC11 CGI-135
Protein Name:	Mitochondrial fission 1 protein (FIS1 homolog) (hFis1) (Tetratricopeptide repeat protein 11) (TPR repeat protein 11)
Other Name:	
Human Swiss Prot No.:	Q9Y3D6
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of human FIS1
Storage and Stability:	-20°C/1 year