

Antibodies Product Datasheet

Product Name:	GAK3 Antibody
Catalog Number:	C42206
Reactivity:	Human, Rat, Mouse
Applications:	WB
Dilutions:	WB 1:1000-2000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human GAK3
Concentration:	1 mg/ml
Gene Name:	
Protein Name:	HERV-K_19q12 provirus ancestral Gag polyprotein (Gag polyprotein) (HERV-K(C19) Gag protein) [Cleaved into: Matrix protein; Capsid protein; Nucleocapsid protein]
Other Name:	
Human Swiss Prot No.:	Q9YNA8
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of human GAK3
Storage and Stability:	-20°C/1 year