

Antibodies Product Datasheet

Product Name:	GAPDH Antibody(Drosophilidae)
Catalog Number:	C42208
Reactivity:	Drosophilidae
Applications:	WB
Dilutions:	WB 1:5000-20000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Purified recombinant Full length GAPDH (Drosophilidae)
Concentration:	1 mg/ml
Gene Name:	
Protein Name:	GAPDH
Other Name:	GAPDH; YM33039
Human Swiss Prot No.:	
Formulation:	Liquid in PBS containing 50% glycerol and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of GAPDH (Drosophilidae)
Storage and Stability:	-20°C/1 year