

Antibodies Product Datasheet

Product Name:	GLGB Antibody
Catalog Number:	C42214
Reactivity:	Human, Mouse
Applications:	WB
Dilutions:	WB 1:1000-2000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human GLGB
Concentration:	1 mg/ml
Gene Name:	GBE1
Protein Name:	1,4-alpha-glucan-branching enzyme (EC 2.4.1.18) (Brancher enzyme) (Glycogen-branching enzyme)
Other Name:	
Human Swiss Prot No.:	Q04446
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of human GLGB
Storage and Stability:	-20°C/1 year