

Antibodies Product Datasheet

Product Name:	GOG6A Antibody
Catalog Number:	C42218
Reactivity:	Human, Rat, Mouse
Applications:	WB
Dilutions:	WB 1:1000-2000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human GOG6A
Concentration:	1 mg/ml
Gene Name:	GOLGA6A GLP GOLGA6
Protein Name:	Golgin subfamily A member 6A (Golgin linked to PML) (Golgin-like protein)
Other Name:	
Human Swiss Prot No.:	Q9NYA3
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of human GOG6A
Storage and Stability:	-20°C/1 year