

Antibodies Product Datasheet

Product Name:	GPDM Antibody
Catalog Number:	C42219
Reactivity:	Human, Mouse, Rat
Applications:	WB
Dilutions:	WB 1:1000-2000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human GPDM
Concentration:	1 mg/ml
Gene Name:	GPD2
Protein Name:	Glycerol-3-phosphate dehydrogenase, mitochondrial (GPD-M) (GPDH-M) (EC 1.1.5.3) (mtGPD)
Other Name:	
Human Swiss Prot No.:	P43304
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of human GPDM
Storage and Stability:	-20°C/1 year