

Antibodies Product Datasheet

Product Name:	GXLT2 Antibody
Catalog Number:	C42223
Reactivity:	Human, Mouse
Applications:	WB
Dilutions:	WB 1:1000-2000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human GXLT2
Concentration:	1 mg/ml
Gene Name:	GXYLT2 GLT8D4
Protein Name:	Glucoside xylosyltransferase 2 (EC 2.4.2.n2) (Glycosyltransferase 8 domain-containing protein 4)
Other Name:	
Human Swiss Prot No.:	A0PJZ3
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of human GXLT2
Storage and Stability:	-20°C/1 year