

Antibodies Product Datasheet

Product Name:	GβL Antibody
Catalog Number:	C42225
Reactivity:	Human, Mouse, Rat
Applications:	IHC-p, WB
Dilutions:	IHC-p1:50-200 ,WB 1:1000-2000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human GβL
Concentration:	1 mg/ml
Gene Name:	MLST8 GBL LST8
Protein Name:	GβL
Other Name:	Target of rapamycin complex subunit LST8 (TORC subunit LST8) (G protein beta subunit-like) (Gable) (Protein GbetaL) (Mammalian lethal with SEC13 protein 8) (mLST8); YM3815
Human Swiss Prot No.:	Q9BVC4
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.87% sodium azide.
Specificity:	This antibody detects endogenous levels of GβL at Human, Mouse, Rat
Storage and Stability:	-20°C/1 year