

Antibodies Product Datasheet

Product Name:	HEM0 Antibody
Catalog Number:	C42230
Reactivity:	Human, Mouse, Rat
Applications:	WB
Dilutions:	WB 1:1000-2000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human HEM0
Concentration:	1 mg/ml
Gene Name:	ALAS2 ALASE ASB
Protein Name:	5-aminolevulinate synthase, erythroid-specific, mitochondrial (ALAS-E) (EC 2.3.1.37) (5-aminolevulinic acid synthase 2) (Delta-ALA synthase 2) (Delta-aminolevulinate synthase 2)
Other Name:	
Human Swiss Prot No.:	P22557
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of human HEM0
Storage and Stability:	-20°C/1 year