

Antibodies Product Datasheet

Product Name:	HM20A Antibody
Catalog Number:	C42234
Reactivity:	Human, Mouse
Applications:	WB
Dilutions:	WB 1:1000-2000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human HM20A
Concentration:	1 mg/ml
Gene Name:	HMG20A HMGX1 HMGXB1
Protein Name:	High mobility group protein 20A (HMG box-containing protein 20A) (HMG domain-containing protein 1) (HMG domain-containing protein HMGX1)
Other Name:	
Human Swiss Prot No.:	Q9NP66
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of human HM20A
Storage and Stability:	-20°C/1 year