

Antibodies Product Datasheet

Product Name:	HOP Antibody
Catalog Number:	C42236
Reactivity:	Human, Mouse, Rat
Applications:	WB
Dilutions:	WB 1:1000-2000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human HOP
Concentration:	1 mg/ml
Gene Name:	HOPX HOD HOP LAGY NECC1 OB1
Protein Name:	Homeodomain-only protein (Lung cancer-associated Y protein) (Not expressed in choriocarcinoma protein 1) (Odd homeobox protein 1)
Other Name:	
Human Swiss Prot No.:	Q9BPY8
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of human HOP
Storage and Stability:	-20°C/1 year