

Antibodies Product Datasheet

Product Name:	HXD4 Antibody
Catalog Number:	C42238
Reactivity:	Human, Mouse
Applications:	WB
Dilutions:	WB 1:1000-2000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human HXD4
Concentration:	1 mg/ml
Gene Name:	HOXD4 HOX4B
Protein Name:	Homeobox protein Hox-D4 (Homeobox protein HHO.C13) (Homeobox protein Hox-4B) (Homeobox protein Hox-5.1)
Other Name:	
Human Swiss Prot No.:	P09016
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of human HXD4
Storage and Stability:	-20°C/1 year