

Antibodies Product Datasheet

Product Name:	ISK1 Antibody
Catalog Number:	C42243
Reactivity:	Human, Rat
Applications:	WB
Dilutions:	WB 1:1000-2000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human ISK1
Concentration:	1 mg/ml
Gene Name:	SPINK1 PSTI
Protein Name:	Pancreatic secretory trypsin inhibitor (Serine protease inhibitor Kazal-type 1) (Tumor-associated trypsin inhibitor) (TATI)
Other Name:	
Human Swiss Prot No.:	P00995
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of human ISK1
Storage and Stability:	-20°C/1 year