

Antibodies Product Datasheet

Product Name:	ISOC2 Antibody
Catalog Number:	C42244
Reactivity:	Human, Rat
Applications:	WB
Dilutions:	WB 1:1000-2000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human ISOC2
Concentration:	1 mg/ml
Gene Name:	ISOC2
Protein Name:	Isochorismatase domain-containing protein 2, mitochondrial
Other Name:	
Human Swiss Prot No.:	Q96AB3
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of human ISOC2
Storage and Stability:	-20°C/1 year