

## Antibodies Product Datasheet

|                               |                                                                                                           |
|-------------------------------|-----------------------------------------------------------------------------------------------------------|
| <b>Product Name:</b>          | K1161 Antibody                                                                                            |
| <b>Catalog Number:</b>        | C42250                                                                                                    |
| <b>Reactivity:</b>            | Human, Mouse                                                                                              |
| <b>Applications:</b>          | WB                                                                                                        |
| <b>Dilutions:</b>             | WB 1:1000-2000                                                                                            |
| <b>Modification:</b>          | Unmodified/Total                                                                                          |
| <b>Source:</b>                | Polyclonal Rabbit                                                                                         |
| <b>Purification:</b>          | The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen. |
| <b>Immunogen:</b>             | Synthesized peptide derived from human K1161                                                              |
| <b>Concentration:</b>         | 1 mg/ml                                                                                                   |
| <b>Gene Name:</b>             | KIAA1161                                                                                                  |
| <b>Protein Name:</b>          | Uncharacterized family 31 glucosidase KIAA1161 (EC 3.2.1.-)                                               |
| <b>Other Name:</b>            |                                                                                                           |
| <b>Human Swiss Prot No.:</b>  | Q6NSJ0                                                                                                    |
| <b>Formulation:</b>           | Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.                                   |
| <b>Specificity:</b>           | This antibody detects endogenous levels of human K1161                                                    |
| <b>Storage and Stability:</b> | -20°C/1 year                                                                                              |