

Antibodies Product Datasheet

Product Name:	KCC2A Antibody
Catalog Number:	C42251
Reactivity:	Human, Mouse, Rat
Applications:	WB
Dilutions:	WB 1:1000-2000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human KCC2A
Concentration:	1 mg/ml
Gene Name:	CAMK2A CAMKA KIAA0968
Protein Name:	Calcium/calmodulin-dependent protein kinase type II subunit alpha (CaM kinase II subunit alpha) (CaMK-II subunit alpha) (EC 2.7.11.17)
Other Name:	
Human Swiss Prot No.:	Q9UQM7
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of human KCC2A
Storage and Stability:	-20°C/1 year