

Antibodies Product Datasheet

Product Name:	KDM4E Antibody
Catalog Number:	C42253
Reactivity:	Human, Rat, Mouse
Applications:	WB
Dilutions:	WB 1:1000-2000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human KDM4E
Concentration:	1 mg/ml
Gene Name:	KDM4E KDM4DL
Protein Name:	Lysine-specific demethylase 4E (EC 1.14.11.-) (KDM4D-like protein) (Lysine-specific demethylase 4D-like)
Other Name:	
Human Swiss Prot No.:	B2RXH2
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of human KDM4E
Storage and Stability:	-20°C/1 year