

Antibodies Product Datasheet

Product Name:	KDM7 Antibody
Catalog Number:	C42254
Reactivity:	Human, Mouse
Applications:	WB
Dilutions:	WB 1:1000-2000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human KDM7
Concentration:	1 mg/ml
Gene Name:	JHDM1D KDM7 KIAA1718
Protein Name:	Lysine-specific demethylase 7 (EC 1.14.11.-) (JmjC domain-containing histone demethylation protein 1D)
Other Name:	
Human Swiss Prot No.:	Q6ZMT4
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of human KDM7
Storage and Stability:	-20°C/1 year