

Antibodies Product Datasheet

Product Name:	KIAS3 Antibody
Catalog Number:	C42256
Reactivity:	Human, Rat, Mouse
Applications:	WB
Dilutions:	WB 1:1000-2000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human KIAS3
Concentration:	1 mg/ml
Gene Name:	KIRREL3-AS3 NCRNA00288 PRR10
Protein Name:	Putative uncharacterized protein KIRREL3-AS3 (KIRREL3 antisense RNA 1) (KIRREL3 antisense gene protein 1) (Proline-rich protein 10)
Other Name:	
Human Swiss Prot No.:	Q8N7Y1
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of human KIAS3
Storage and Stability:	-20°C/1 year