

Antibodies Product Datasheet

Product Name:	KIFA3 Antibody
Catalog Number:	C42258
Reactivity:	Human, Mouse
Applications:	WB
Dilutions:	WB 1:1000-2000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human KIFA3
Concentration:	1 mg/ml
Gene Name:	KIFAP3 KIF3AP SMAP
Protein Name:	Kinesin-associated protein 3 (KAP-3) (KAP3) (Smg GDS-associated protein)
Other Name:	
Human Swiss Prot No.:	Q92845
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of human KIFA3
Storage and Stability:	-20°C/1 year