

Antibodies Product Datasheet

Product Name:	Kininogen 1 Rabbit Antibody
Catalog Number:	C42259
Reactivity:	Human, Mouse, Rat
Applications:	WB, IHC
Dilutions:	IHC-p 1:50-200, WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Recombinant Protein derived from human Kininogen 1
Concentration:	1 mg/ml
Gene Name:	KNG1 BDK KNG
Protein Name:	Kininogen 1
Other Name:	Kininogen-1 (Alpha-2-thiol proteinase inhibitor) (Fitzgerald factor) (High molecular weight kininogen) (HMWK) (Williams-Fitzgerald-Flaujeac factor) [Cleaved into: Kininogen-1 heavy chain; T-kinin (Ile-Ser-Bradykinin); Bradykinin (Kallidin I); Lysyl-bradyk
Human Swiss Prot No.:	P01042
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of human Kininogen 1
Storage and Stability:	-20°C/1 year