

Antibodies Product Datasheet

Product Name:	KLH31 Antibody
Catalog Number:	C42261
Reactivity:	Human, Mouse
Applications:	WB
Dilutions:	WB 1:1000-2000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human KLH31
Concentration:	1 mg/ml
Gene Name:	KLHL31 BKLHD6 KBTBD1 KLHL
Protein Name:	Kelch-like protein 31 (BTB and kelch domain-containing protein 6) (Kelch repeat and BTB domain-containing protein 1) (Kelch-like protein KLHL)
Other Name:	
Human Swiss Prot No.:	Q9H511
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of human KLH31
Storage and Stability:	-20°C/1 year