

Antibodies Product Datasheet

Product Name:	LAP4A Antibody
Catalog Number:	C42268
Reactivity:	Human, Mouse, Rat
Applications:	WB
Dilutions:	WB 1:1000-2000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human LAP4A
Concentration:	1 mg/ml
Gene Name:	LAPTM4A KIAA0108 LAPTM4 MBNT MTRP UNQ1846/PRO3574
Protein Name:	Lysosomal-associated transmembrane protein 4A (Golgi 4-transmembrane-spanning transporter MTP)
Other Name:	
Human Swiss Prot No.:	Q15012
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of human LAP4A
Storage and Stability:	-20°C/1 year