

Antibodies Product Datasheet

Product Name:	LMAN2 Antibody
Catalog Number:	C42275
Reactivity:	Human, Mouse
Applications:	WB
Dilutions:	WB 1:1000-2000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human LMAN2
Concentration:	1 mg/ml
Gene Name:	LMAN2 C5orf8
Protein Name:	Vesicular integral-membrane protein VIP36 (Glycoprotein GP36b) (Lectin mannose-binding 2) (Vesicular integral-membrane protein 36) (VIP36)
Other Name:	
Human Swiss Prot No.:	Q12907
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of human LMAN2
Storage and Stability:	-20°C/1 year