

Antibodies Product Datasheet

Product Name:	LMBL3 Antibody
Catalog Number:	C42276
Reactivity:	Human, Mouse
Applications:	WB
Dilutions:	WB 1:1000-2000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human LMBL3
Concentration:	1 mg/ml
Gene Name:	L3MBTL3 KIAA1798 MBT1
Protein Name:	Lethal(3)malignant brain tumor-like protein 3 (H-I(3)mbt-like protein 3) (L(3)mbt-like protein 3) (MBT-1)
Other Name:	
Human Swiss Prot No.:	Q96JM7
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of human LMBL3
Storage and Stability:	-20°C/1 year