

Antibodies Product Datasheet

Product Name:	MA1A2 Antibody
Catalog Number:	C42283
Reactivity:	Human, Mouse
Applications:	WB
Dilutions:	WB 1:1000-2000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human MA1A2
Concentration:	1 mg/ml
Gene Name:	MAN1A2 MAN1B
Protein Name:	Mannosyl-oligosaccharide 1,2-alpha-mannosidase IB (EC 3.2.1.113) (Mannosidase alpha class 1A member 2) (Processing alpha-1,2-mannosidase IB) (Alpha-1,2-mannosidase IB)
Other Name:	
Human Swiss Prot No.:	O60476
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of human MA1A2
Storage and Stability:	-20°C/1 year