

Antibodies Product Datasheet

Product Name:	MMSA Antibody
Catalog Number:	C42299
Reactivity:	Human, Mouse, Rat
Applications:	WB
Dilutions:	WB 1:1000-2000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human MMSA
Concentration:	1 mg/ml
Gene Name:	ALDH6A1 MMSDH
Protein Name:	Methylmalonate-semialdehyde dehydrogenase [acylating], mitochondrial (MMSDH) (Malonate-semialdehyde dehydrogenase [acylating]) (EC 1.2.1.18) (EC 1.2.1.27) (Aldehyde dehydrogenase family 6 member A1)
Other Name:	
Human Swiss Prot No.:	Q02252
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of human MMSA
Storage and Stability:	-20°C/1 year