

## Antibodies Product Datasheet

|                               |                                                                                                                                                                                                                                                                |
|-------------------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Product Name:</b>          | NDST3 Antibody                                                                                                                                                                                                                                                 |
| <b>Catalog Number:</b>        | C42322                                                                                                                                                                                                                                                         |
| <b>Reactivity:</b>            | Human, Mouse                                                                                                                                                                                                                                                   |
| <b>Applications:</b>          | WB                                                                                                                                                                                                                                                             |
| <b>Dilutions:</b>             | WB 1:1000-2000                                                                                                                                                                                                                                                 |
| <b>Modification:</b>          | Unmodified/Total                                                                                                                                                                                                                                               |
| <b>Source:</b>                | Polyclonal Rabbit                                                                                                                                                                                                                                              |
| <b>Purification:</b>          | The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.                                                                                                                                                      |
| <b>Immunogen:</b>             | Synthesized peptide derived from human NDST3                                                                                                                                                                                                                   |
| <b>Concentration:</b>         | 1 mg/ml                                                                                                                                                                                                                                                        |
| <b>Gene Name:</b>             | NDST3 HSST3 UNQ2544/PRO4998                                                                                                                                                                                                                                    |
| <b>Protein Name:</b>          | Bifunctional heparan sulfate N-deacetylase/N-sulfotransferase 3 (EC 2.8.2.8) (Glucosaminyl N-deacetylase/N-sulfotransferase 3) (NDST-3) (hNDST-3) (N-heparan sulfate sulfotransferase 3) (N-HSST 3) [Includes: Heparan sulfate N-deacetylase 3 (EC 3.-.-.); He |
| <b>Other Name:</b>            |                                                                                                                                                                                                                                                                |
| <b>Human Swiss Prot No.:</b>  | O95803                                                                                                                                                                                                                                                         |
| <b>Formulation:</b>           | Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.                                                                                                                                                                                        |
| <b>Specificity:</b>           | This antibody detects endogenous levels of human NDST3                                                                                                                                                                                                         |
| <b>Storage and Stability:</b> | -20°C/1 year                                                                                                                                                                                                                                                   |