

Antibodies Product Datasheet

Product Name:	NEUR1 Antibody
Catalog Number:	C42323
Reactivity:	Human, Mouse, Rat
Applications:	WB
Dilutions:	WB 1:1000-2000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human NEUR1
Concentration:	1 mg/ml
Gene Name:	NEU1 NANH
Protein Name:	Sialidase-1 (EC 3.2.1.18) (Acetylneuraminyl hydrolase) (G9 sialidase) (Lysosomal sialidase) (N-acetyl-alpha-neuraminidase 1)
Other Name:	
Human Swiss Prot No.:	Q99519
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of human NEUR1
Storage and Stability:	-20°C/1 year