

Antibodies Product Datasheet

Product Name:	NMNA2 Antibody
Catalog Number:	C42326
Reactivity:	Human, Mouse, Rat
Applications:	WB
Dilutions:	WB 1:1000-2000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human NMNA2
Concentration:	1 mg/ml
Gene Name:	NMNAT2 C1orf15 KIAA0479
Protein Name:	Nicotinamide mononucleotide adenylyltransferase 2 (NMN adenylyltransferase 2) (EC 2.7.7.1) (Nicotinate-nucleotide adenylyltransferase 2) (NaMN adenylyltransferase 2) (EC 2.7.7.18)
Other Name:	
Human Swiss Prot No.:	Q9BZQ4
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of human NMNA2
Storage and Stability:	-20°C/1 year