
Antibodies Product Datasheet

Product Name:	1A26 Antibody
Catalog Number:	C42379
Reactivity:	Human, Rat, Mouse
Applications:	WB
Dilutions:	WB 1: 500-2000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human 1A26
Concentration:	1 mg/ml
Gene Name:	HLA-A HLAA
Protein Name:	1A26
Other Name:	
Human Swiss Prot No.:	P30450
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of 1A26 at Human
Storage and Stability:	-20°C/1 year