

Antibodies Product Datasheet

Product Name:	2A5G Antibody
Catalog Number:	C42388
Reactivity:	Human, Mouse
Applications:	WB, ELISA
Dilutions:	WB 1:500-2000 ELISA 1:5000-20000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human protein . at AA range: 360-440
Concentration:	1 mg/ml
Gene Name:	PPP2R5C KIAA0044
Protein Name:	Serine/threonine-protein phosphatase 2A 56 kDa regulatory subunit gamma isoform (PP2A B subunit isoform B'-gamma) (PP2A B subunit isoform B56-gamma) (PP2A B subunit isoform PR61-gamma) (PP2A B subunit isoform R5-gamma) (Renal carcinoma antigen NY-REN-29)
Other Name:	
Human Swiss Prot No.:	Q13362
Formulation:	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Specificity:	2A5G Polyclonal Antibody detects endogenous levels of protein.
Storage and Stability:	-20°C/1 year