

Antibodies Product Datasheet

Product Name:	2AAA Antibody
Catalog Number:	C42389
Reactivity:	Human, Mouse
Applications:	WB, ELISA
Dilutions:	WB 1:500-2000 ELISA 1:5000-20000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human protein . at AA range: 500-580
Concentration:	1 mg/ml
Gene Name:	PPP2R1A
Protein Name:	Serine/threonine-protein phosphatase 2A 65 kDa regulatory subunit A alpha isoform (Medium tumor antigen-associated 61 kDa protein) (PP2A subunit A isoform PR65-alpha) (PP2A subunit A isoform R1-alpha)
Other Name:	
Human Swiss Prot No.:	P30153
Formulation:	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Specificity:	2AAA Polyclonal Antibody detects endogenous levels of protein.
Storage and Stability:	-20°C/1 year