
Antibodies Product Datasheet

| | |
|-------------------------------|---|
| Product Name: | 2B1F Antibody |
| Catalog Number: | C42393 |
| Reactivity: | Human, Rat, Mouse |
| Applications: | WB, ELISA |
| Dilutions: | WB 1:500-2000 ELISA 1:5000-20000 |
| Modification: | Unmodified/Total |
| Source: | Polyclonal Rabbit |
| Purification: | The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen. |
| Immunogen: | Synthesized peptide derived from part region of human protein |
| Concentration: | 1 mg/ml |
| Gene Name: | HLA-DRB1 HLA-DRB2 |
| Protein Name: | HLA class II histocompatibility antigen, DRB1-15 beta chain (DW2.2/DR2.2) (MHC class II antigen DRB1*15) |
| Other Name: | |
| Human Swiss Prot No.: | P01911 |
| Formulation: | Liquid in PBS containing 50% glycerol, and 0.02% sodium azide. |
| Specificity: | 2B1F Polyclonal Antibody detects endogenous levels of protein. |
| Storage and Stability: | -20°C/1 year |