

Antibodies Product Datasheet

Product Name:	3BHS1 Antibody
Catalog Number:	C42394
Reactivity:	Human, Rat, Mouse
Applications:	WB, ELISA
Dilutions:	WB 1:500-2000 ELISA 1:5000-20000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human protein . at AA range: 30-110
Concentration:	1 mg/ml
Gene Name:	HSD3B1 3BH HSDB3A
Protein Name:	3 beta-hydroxysteroid dehydrogenase/Delta 5-->4-isomerase type 1 (3 beta-hydroxysteroid dehydrogenase/Delta 5-->4-isomerase type I) (3-beta-HSD I) (Trophoblast antigen FDO161G) [Includes: 3-beta-hydroxy-Delta(5)-steroid dehydrogenase (EC 1.1.1.145) (3-bet
Other Name:	
Human Swiss Prot No.:	P14060
Formulation:	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Specificity:	3BHS1 Polyclonal Antibody detects endogenous levels of protein.
Storage and Stability:	-20°C/1 year