

## Antibodies Product Datasheet

<b>Product Name:</b>	3BHS2 Antibody
<b>Catalog Number:</b>	C42395
<b>Reactivity:</b>	Human, Rat, Mouse
<b>Applications:</b>	WB, ELISA
<b>Dilutions:</b>	WB 1:500-2000 ELISA 1:5000-20000
<b>Modification:</b>	Unmodified/Total
<b>Source:</b>	Polyclonal Rabbit
<b>Purification:</b>	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
<b>Immunogen:</b>	Synthesized peptide derived from human protein . at AA range: 180-260
<b>Concentration:</b>	1 mg/ml
<b>Gene Name:</b>	HSD3B2 HSDB3B
<b>Protein Name:</b>	3 beta-hydroxysteroid dehydrogenase/Delta 5-->4-isomerase type 2 (3 beta-hydroxysteroid dehydrogenase/Delta 5-->4-isomerase type II) (3-beta-HSD II) (3-beta-HSD adrenal and gonadal type) [Includes: 3-beta-hydroxy-Delta(5)-steroid dehydrogenase (EC 1.1.1.1
<b>Other Name:</b>	
<b>Human Swiss Prot No.:</b>	P26439
<b>Formulation:</b>	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
<b>Specificity:</b>	3BHS2 Polyclonal Antibody detects endogenous levels of protein.
<b>Storage and Stability:</b>	-20°C/1 year