

Antibodies Product Datasheet

Product Name:	5NTC Antibody
Catalog Number:	C42404
Reactivity:	Human, Mouse
Applications:	WB, ELISA
Dilutions:	WB 1:500-2000 ELISA 1:5000-20000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from part region of human protein
Concentration:	1 mg/ml
Gene Name:	NT5C2 NT5B NT5CP PNT5
Protein Name:	Cytosolic purine 5'-nucleotidase (EC 3.1.3.5) (Cytosolic 5'-nucleotidase II)
Other Name:	
Human Swiss Prot No.:	P49902
Formulation:	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Specificity:	5NTC Polyclonal Antibody detects endogenous levels of protein.
Storage and Stability:	-20°C/1 year