

Antibodies Product Datasheet

Product Name:	7B2 Antibody
Catalog Number:	C42405
Reactivity:	Human, Mouse, Rat
Applications:	WB, ELISA
Dilutions:	WB 1:500-2000 ELISA 1:5000-20000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human protein . at AA range: 90-170
Concentration:	1 mg/ml
Gene Name:	SCG5 SGNE1
Protein Name:	Neuroendocrine protein 7B2 (Pituitary polypeptide) (Secretogranin V) (Secretogranin-5) (Secretory granule endocrine protein I) [Cleaved into: N-terminal peptide; C-terminal peptide]
Other Name:	
Human Swiss Prot No.:	P05408
Formulation:	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Specificity:	7B2 Polyclonal Antibody detects endogenous levels of protein.
Storage and Stability:	-20°C/1 year