

Antibodies Product Datasheet

Product Name:	A4GAT Antibody
Catalog Number:	C42412
Reactivity:	Human, Rat, Mouse
Applications:	WB, ELISA
Dilutions:	WB 1:500-2000 ELISA 1:5000-20000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from part region of human protein
Concentration:	1 mg/ml
Gene Name:	A4GALT A14GALT A4GALT1
Protein Name:	Lactosylceramide 4-alpha-galactosyltransferase (EC 2.4.1.228) (Alpha-1,4-N-acetylglucosaminyltransferase) (Alpha-1,4-galactosyltransferase) (Alpha4Gal-T1) (CD77 synthase) (Globotriaosylceramide synthase) (Gb3 synthase) (P1/Pk synthase) (UDP-galactose:beta
Other Name:	
Human Swiss Prot No.:	Q9NPC4
Formulation:	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Specificity:	A4GAT Polyclonal Antibody detects endogenous levels of protein.
Storage and Stability:	-20°C/1 year