

Antibodies Product Datasheet

Product Name:	ABC3C Antibody
Catalog Number:	C42419
Reactivity:	Human, Rat, Mouse
Applications:	WB, ELISA
Dilutions:	WB 1:500-2000 ELISA 1:5000-20000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human protein . at AA range: 10-90
Concentration:	1 mg/ml
Gene Name:	APOBEC3C APOBEC1L PBI
Protein Name:	Probable DNA dC->dU-editing enzyme APOBEC-3C (EC 3.5.4.-) (APOBEC1-like) (Phorbolin I)
Other Name:	
Human Swiss Prot No.:	Q9NRW3
Formulation:	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Specificity:	ABC3C Polyclonal Antibody detects endogenous levels of protein.
Storage and Stability:	-20°C/1 year