

Antibodies Product Datasheet

Product Name:	ABCAC Antibody
Catalog Number:	C42426
Reactivity:	Human, Rat, Mouse
Applications:	IHC-p, IF(paraffin section)
Dilutions:	IHC-p 1:50-300
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human protein . at AA range: 2170-2250
Concentration:	1 mg/ml
Gene Name:	ABCA12 ABC12
Protein Name:	ATP-binding cassette sub-family A member 12 (ATP-binding cassette transporter 12) (ATP-binding cassette 12)
Other Name:	
Human Swiss Prot No.:	Q86UK0
Formulation:	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Specificity:	ABCAC Polyclonal Antibody detects endogenous levels of protein.
Storage and Stability:	-20°C/1 year