

Antibodies Product Datasheet

Product Name:	ABCD2 Antibody
Catalog Number:	C42433
Reactivity:	Human, Rat, Mouse
Applications:	WB, ELISA
Dilutions:	WB 1:500-2000 ELISA 1:5000-20000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human protein . at AA range: 310-390
Concentration:	1 mg/ml
Gene Name:	ABCD2 ALD1 ALDL1 ALDR ALDRP
Protein Name:	ATP-binding cassette sub-family D member 2 (Adrenoleukodystrophy-like 1) (Adrenoleukodystrophy-related protein) (hALDR)
Other Name:	
Human Swiss Prot No.:	Q9UBJ2
Formulation:	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Specificity:	ABCD2 Polyclonal Antibody detects endogenous levels of protein.
Storage and Stability:	-20°C/1 year