

Antibodies Product Datasheet

Product Name:	ABCD3 Antibody
Catalog Number:	C42434
Reactivity:	Human, Mouse, Rat
Applications:	WB, ELISA
Dilutions:	WB 1:500-2000 ELISA 1:5000-20000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human protein . at AA range: 40-120
Concentration:	1 mg/ml
Gene Name:	ABCD3 PMP70 PXMP1
Protein Name:	ATP-binding cassette sub-family D member 3 (70 kDa peroxisomal membrane protein) (PMP70)
Other Name:	
Human Swiss Prot No.:	P28288
Formulation:	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Specificity:	ABCD3 Polyclonal Antibody detects endogenous levels of protein.
Storage and Stability:	-20°C/1 year