

Antibodies Product Datasheet

Product Name:	ACATN Antibody
Catalog Number:	C42451
Reactivity:	Human, Rat, Mouse
Applications:	WB, ELISA
Dilutions:	WB 1:500-2000 ELISA 1:5000-20000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human protein . at AA range: 110-190
Concentration:	1 mg/ml
Gene Name:	SLC33A1 ACATN AT1
Protein Name:	Acetyl-coenzyme A transporter 1 (AT-1) (Acetyl-CoA transporter 1) (Solute carrier family 33 member 1)
Other Name:	
Human Swiss Prot No.:	O00400
Formulation:	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Specificity:	ACATN Polyclonal Antibody detects endogenous levels of protein.
Storage and Stability:	-20°C/1 year