

Antibodies Product Datasheet

Product Name:	ACL6B Antibody
Catalog Number:	C42462
Reactivity:	Human, Mouse
Applications:	WB, ELISA
Dilutions:	WB 1:500-2000 ELISA 1:5000-20000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human protein . at AA range: 1-80
Concentration:	1 mg/ml
Gene Name:	ACTL6B ACTL6 BAF53B
Protein Name:	Actin-like protein 6B (53 kDa BRG1-associated factor B) (Actin-related protein Baf53b) (ArpNalpha) (BRG1-associated factor 53B) (BAF53B)
Other Name:	
Human Swiss Prot No.:	O94805
Formulation:	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Specificity:	ACL6B Polyclonal Antibody detects endogenous levels of protein.
Storage and Stability:	-20°C/1 year