

Antibodies Product Datasheet

Product Name:	ACO11 Antibody
Catalog Number:	C42464
Reactivity:	Human, Mouse
Applications:	WB, ELISA
Dilutions:	WB 1:500-2000 ELISA 1:5000-20000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human protein . at AA range: 290-370
Concentration:	1 mg/ml
Gene Name:	ACOT11 BFIT KIAA0707 THEA
Protein Name:	Acyl-coenzyme A thioesterase 11 (Acyl-CoA thioesterase 11) (EC 3.1.2.-) (Acyl-CoA thioester hydrolase 11) (Adipose-associated thioesterase) (Brown fat-inducible thioesterase) (BFIT)
Other Name:	
Human Swiss Prot No.:	Q8WXL4
Formulation:	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Specificity:	ACO11 Polyclonal Antibody detects endogenous levels of protein.
Storage and Stability:	-20°C/1 year