

Antibodies Product Datasheet

Product Name:	ACO13 Antibody
Catalog Number:	C42465
Reactivity:	Human, Mouse
Applications:	WB
Dilutions:	WB 1:1000-2000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human ACO13
Concentration:	1 mg/ml
Gene Name:	ACOT13 THEM2 HT012 PNAS-27
Protein Name:	Acyl-coenzyme A thioesterase 13 (Acyl-CoA thioesterase 13) (EC 3.1.2.-) (Thioesterase superfamily member 2)
Other Name:	
Human Swiss Prot No.:	Q9NPJ3
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of human ACO13
Storage and Stability:	-20°C/1 year