

Antibodies Product Datasheet

Product Name:	ACOX2 Antibody
Catalog Number:	C42468
Reactivity:	Human, Mouse, Rat
Applications:	WB, ELISA
Dilutions:	WB 1:500-2000 ELISA 1:5000-20000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human protein . at AA range: 270-350
Concentration:	1 mg/ml
Gene Name:	ACOX2
Protein Name:	Peroxisomal acyl-coenzyme A oxidase 2 (EC 1.17.99.3) (3- α ,7- α ,12- α -trihydroxy-5- β -cholestanoyl-CoA 24-hydroxylase) (3- α ,7- α ,12- α -trihydroxy-5- β -cholestanoyl-CoA oxidase) (Trihydroxycoprostanoyl-CoA oxidase) (THCA-CoA oxidase)
Other Name:	
Human Swiss Prot No.:	Q99424
Formulation:	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Specificity:	ACOX2 Polyclonal Antibody detects endogenous levels of protein.
Storage and Stability:	-20°C/1 year